\land CONFIGURE

Set up Operating Modes

Use these commands to set up the operating configuration of your Model 600A/620A Environmental Chamber System Controller.

C O M M A N D	FUNCTION REPLY
Use these commands to first select the 600/	7620 you wish to address and place the selected unit under COMPUTER control:
ADR [cr]	read the address of the selected unit nn [lf] Use this QUERY command if your installation has two or more 600/620s
#nn [cr]	select the 600/620 ADDRESS [no reply]
IFS [cr]	request INTERFACE STATUS LCL [If], RMT [If], or RSR [If] Control is either from the Front Panel ("LCL") or from your computer ("RMT" or "RSR")
RMT SELECT [cr]	change to REMOTE mode (computer control)ok [If] To transfer control back to the 600/620 Front Panel, use the command LCL [cr]
RSR SELECT [cr]	change to EXTERNAL RAMP SOAK (computer control) ok [If]
2 Use the following commands to READ inform	IATION ABOUT THE SELECTED UNIT:
CRM [cr]	read CONTROLLER MODE
FRE [cr]	read FREE MEMORY
MOD [cr]	read PROGRAMMER MODE
PFA [ɑr]	read POWER FAIL RECOVERY ACTION
I8A [cr]	read #8 INPUT JUMP CONDITION
STA [cr]	read STATUS OF PROGRAMMER
VER [cr]	read FIRMWARE VERSION
DAT [cr]	read DATE of Firmware Revision
SYN [cr]*	read STATUS OF SYNCHRONIZER MODE x [If] Response "x" is one of the following: "1" = ON, "0" = OFF

continued on next page

* Note: Sync must be OFF unless 600A/620A is being used as part of a multi-unit synchronized system!

▲ CONFIGURE (continued)

Set up Operating Modes

Use these commands to set up the operating configuration of your Model 600A/620A Environmental Chamber System Controller.

$\mathsf{C} \ \mathsf{O} \ \mathsf{M} \ \mathsf{M} \ \mathsf{A} \ \mathsf{N} \ \mathsf{D}$

FUNCTION

REPLY

Use these commands	TO CONFIGURE THE SE	LECTED UNIT: (NOTE: SEND COMMAND W/O DATA TO READ CURRENT VALUE)	
ALT nnnn [cr]		select SITE ELEVATION (altitude, in feet)	ok [lf]
CH1 M [cr]		select MANUAL MODE for Channel 1	ok [lf]
		Use "CH2 M" to select Manual Mode for Channel 2	
CH1 P [cr]		select PROGRAMMER MODE for Channel 1 Use "CH2 P" to select Programmer Mode for Channel 2	ok [lf]
FTH nnn [cr]		set the FastTRAC HI TEMP setpoint limit	ok [lf]
FTL nnn [cr]		set the FastTRAC LO TEMP setpoint limit	ok [lf]
HL1 nnn [cr]		set HI PROCESS LIMIT for Channel 1 Use "HL2 nnn" to set Hi Process Limit for Channel 2	ok [lf]
HS1 nnn [cr]		set HI SPAN LIMIT for Channel 1 Use "HS2 nnn" to set Hi Span Limit for Channel 2	ok [lf]
I8A N [cr]		set #8 Input Jump Condition to NO ACTION This is the Normal (default) condition	ok [lf]
18A R [cr]		set Input #8 Jump Condition to RUN Follow this with command defining Program to be run (see next)	ok [lf]
18P nn [cr]		set PROGRAM NUMBER to run if Jump is executed	ok [lf]
18S nnn [cr]		set STARTING STEP in Selected Program if Jump is executed	ok [lf]
LL1 nnn [cr]		set LOWER PROCESS LIMIT for Channel 1	ok [lf]
LP1 nn [cr]		set LOW PASS FILTER for Channel 1	ok [lf] 1 <i>2</i>
LS1 nnn [cr]		set LOW SPAN LIMIT for Channel 1 Use "LS2 nnn" to set Low Span Limit for Channel 2	ok [lf]
PF1 nnn [cr]		set the POWER FAIL Process Change Limit for Channel 1	ok [lf]
PFP nn [cr]		set POWER FAILURE PROGRAM Number	ok [lf]
PFS nnn [cr]		set STARTING STEP in Selected Power Fail Program	ok [lf]
UN1 nnn [cr]		set UNITS for Channel 1	ok [lf]
CVAL - [-]*		"nnn" is one of the following:F-Degrees-F;C-Degrees-C; RH-RH/WB; LIN-Linear Use "UN2 nnn" command to set Units for Channel 2	
SYN n [cr]^		Set SYNGHKUNIZEK MUDE Use "SYN 1" to turn ON the Synchronizaer mode Use "SYN 0" to turn OFF the Synchronizaer mode	ok [It]

* Note: Sync must be OFF unless 600A/620A is being used as part of a multi-unit synchronized system!

EDIT

Use these commands to 600A/620A Environm	o Create, Modify or Review programs stored in your nental Chamber System Controller.
C O M M A N D	FUNCTION REPLY
1 Use these commands to first select t	HE 600/620 YOU WISH TO ADDRESS AND PLACE THE SELECTED UNIT UNDER COMPUTER CONTROL:
ADR [cr]	
#nn [cr] IFS [cr]	select the 600/620 ADDRESS
RMT SELECT [cr]	change to REMOTE mode (computer control) ok [If] To transfer control back to the 600/620 Front Panel, use the command LCL [cr]
2 If you are generating a NEW program	<i>n</i> , start with this command and then go to Item 4 below:
FRS nn [cr]	insert the FIRST STEP of a new program ok [If] "nn" is the number of the new program
¹³ If you wish to Edit an EXISTING pro	igram, start by selecting the desired PROGRAM and STEP then go to Item 4 below:
PRN nn [cr]	
STN nnn [cr]	choose the STEP number in the selected programok [If] To READ the current STEP number, just enter STN [cr] nnn [If]
Use these commands to Edit the TYP	E of the selected step, then go to Item 5 below:
TYP [cr]	
TYP S [cr] TYP E [cr]	define selected step as STD (Standard = RAMP/SOAK) ok [If] define selected step as EOP (End of Program) ok [if] If your program calls for multiple cycles, follow this command with the LOP command (see below)
TYP L [cr]	define selected step as LOOP
LOP cccc sss [cr]	LOOP from the selected step, a total of cccc cycles back to step sss ok [If]
TYP G [cr]	define selected step as GTO (Go To - Unconditional) ok [If] Immediately follow this command with the GTO command (see next)
GTO pp sss [cr]	directs the controller to JUMP to Program "pp", Step "sss" ok [If]
TYP I [cr]	define selected step as GTI (Go To -IF) ok [If] Immediately follow this command with the GTI command (see next)
GTI pp sss x [cr]	directs the controller to JUMP to Program "pp", Step "sss" IF "x" is true ok [If] "x" may be any event ("1" thru" 8" or "A "thru" H") or "0" for Input #4
TYP P [cr]	define selected step as PAUSE

continued on next page \square

 I E E E 4 8 8
 O P T I O N S

 Refer to the back cover of this guide for a listing of additional IEEE488 instructions.

▲ EDIT (continued)

Create and/or Modify or Review Programs

Use these commands to Create, Modify or Review programs stored in your 600A/620A Environmental Chamber System Controller.

$\mathsf{C} \hspace{0.1cm} \mathsf{O} \hspace{0.1cm} \mathsf{M} \hspace{0.1cm} \mathsf{M} \hspace{0.1cm} \mathsf{A} \hspace{0.1cm} \mathsf{N} \hspace{0.1cm} \mathsf{D}$

FUNCTION

REPLY

	Use these commands t	D COMPLETE THE DETA	ILS for the selected step. (note: send command w/o data to read c	URRENT VALUE)
	AL1 nn [cr]		set DEVIATION ALARM for Channel 1 (use "AL2" for Ch2)	ok [lf]
	DEP [cr]		DELETE the selected PROGRAM	ok [lf]
	DES [cr]		DELETE the selected STEP	ok [lf]
	CL1 [cr]		CLEAR the Channel 1 SETPOINT for the selected step (use "CL2" for Ch2)	ok [lf]
	EVC [cr]		CLEAR all the EVENTS for the selected step	ok [lf]
	EVD [cr]		set ALL the EVENTS for the selected step to DON'T CARE state	ok [lf]
	EVF n n [cr]		set selected EVENTS to OFF (e.g. EVF 1 3 A H [cr])	ok [lf]
	EVN n n [cr]		set selected EVENTS to ON (e.g. EVN 2 4 B G [cr])	ok [lf]
	EVS 1•3••6••A••D•	X•H [cr]	EVENT STATUS (read or store image)	ok [lf]
	EVX n n [cr]		set selected EVENTS to DON'T CARE (e.g. EVX 2 5 C [cr])	ok [lf]
	GS1 nn [cr]		set GUARANTEED SOAK for ch1 (use "GS2" for ch2)	ok [lf]
	SS1 nnn.n [cr]	establish SETPOINT for Channel 1 (use "SS2" for Ch2)	ok [lf]
	STM hh mm ss	[cr]	set TIME for step	ok [lf]
			Enter hours as "00" if time less than 1 hr; You may also enter this command as STM h	nh:mm:ss [cr]
	TBD nn [cr]		set thermal boost decrease	ok [lf]
	TBI nn [cr]		set thermal boost increase	ok [lf]
6	If you are Creating a I	NEW program or ADDIN	IG steps to an existing program, use this command then return to iti	em 4 above:
	INS [cr]		INSERT a BLANK STEP immediately following the selected step	ok [lf]
			All subsequent commands will apply to the NEW step	
7	Use this command to I	RETURN THE SELECTED 60	00/620 Controller to FRONT PANEL operation:	
	LCL [cr]		place selected unit under LOCAL control (front panel)	ok [lf]

I E E E A 8 8 O P T I O N S Refer to the back cover of this guide for a listing of additional IEEE488 instructions. **MONITOR**

Read Current Values & Status while the Program is Running

Use these commands to monitor the current information in your Model 600A/620A Environmental Chamber System Controller.

C O M M A N D	F U N C T I O N R E P L Y
Use these commands to first select the 600/	/620 you wish to address:
ADR [cr]	read the address of the selected unit
#nn [cr]	select the 600/620 ADDRESS
2 Use these commands to request SET UP infor	rmation from the selected 600/620 Controller:
IFS [cr]	request INTERFACE STATUS LCL [IF] or RMT [IF]
CH1 [cr]	Response will be "LCL" if selected unit is under front panel control; "RMT" if under computer control request Channel 1 CONTROLLER MODE
UN1 [α]	read UNITS of Channel 1 Configuration
3 Use these commands to request CONTROL LC	OOP information from the selected 600/620:
PR1 [cr]	read PROCESS VALUE for Channel 1
PW1 [cr]	read POWER OUTPUT for Channel 1 xx xx [If]
	Response is as percentage of maximum
	Use command "PW2" to request Channel 2 Power Output
3P1 [0]	Use command "SP2" to request Channel 2 Setpoint Value
4 Use these commands to request ALARM and E	VENT INFORMATION FROM THE SELECTED 600/620:
AEV [cr]	EVENT STATUS Request 1•••5•7•.A•••E••H [If]
	Events that are "ON" are noted by identifying character; events that are "OFF" are noted by "•"
ALIVI [0]	ALARM STATUS Request
CPI [cr]	read LATCHED inputs (on TB-5) and RESET
RDI [cr]	read STATUS OF INPUTS at TB-6

continued on next page \Box

I E E E 4 8 8O P T I O N SRefer to the back cover of this guide for a listing of additional IEEE488 instructions.

▲ MONITOR (continued)

COMMAND

REPLY

Use these commands to monitor the current information in your Model 600A/620A Environmental Chamber System Controller.

FUNCTION

5 Use these commands to request PROGRAMMER information from the selected 600/620: CPN [cr] read Current PROGRAM Number nn [lf] CSN [cr] read Current STEP Number nnn [lf] MOD [cr] identify MODE OF PROGRAMMER xxx [lf] Response "xxx" is "RUN" (Run Mode) or "STP" (Stop Mode) nnnn [lf] CTG [cr] read CYCLES TO GO nnnn [lf] Response to this command will be "----" if EOP has not been reached nnnn [lf] ITM [cr] read TIME REMAINING in Current Step. hh:mm:ss [lf]

IEEE488 OPTIONS *Refer to the back cover of this guide for a listing of additional IEEE488 instructions.* ▲ PID TUNING

Change Tuning Parameters in Real Time

Use these commands to set or modify the PID tuning parameters of your Model 600A/620A Environmental Chamber System Controller.

C O M M A N D	FUNCTION	REPLY
Use these commands to first select t	THE 600/620 YOU WISH TO ADDRESS AND PLACE THE SEL	ECTED UNIT UNDER COMPUTER CONTROL:
ADR [cr]	read the ADDRESS of the selected unit Use this QUERY command if your installation has two or	nn [lf] more 600/620s
#nn [cr]	select the 600/620 ADDRESS	[no reply]
IFS [cr]		LCL [If] or RMT [If]
RMT SELECT [cr]	change to REMOTE mode (computer control) To transfer control back to the 600/620 Front Panel, us	se the command LCL [cr]
Use these commands to ADDRESS (or read) a group of PID parameters:	
ΡΙΟ nnn [cr]	ADDRESS PID parameters "nnn" is a three-character code as follows: 1 St Character: the CHANNEL NUMBER (either "1" or "2 2 nd Character: HEATING or COOLING action (either "H" 3 rd Character: UNSTAGED or STAGED GROUP (either "U as illustrated in the following example: DID 1444 [cr] Scleet DID (bearned 1, Macting Action UK	or "C") or "C") " or "S")
PID [cr]	READ the active PID parameter selection "nnn" is a three character code as described above	nnn [lf]
Use these commands to SET the PID	PARAMETERS WITHIN THE SELECTED GROUP: (note: for a qu	ery, simply omit the "n" from the command)
CTM n [cr]	set CYCLE TIME	ok [lf]
PBN n [cr]	set PROPORTIONAL BAND	ok [lf]
PGN n [cr]	set PROPORTIONAL GAIN	ok [lf]
INT n [cr]	set Integral (RESET) action	ok [lf]
RWI n [cr]	set RESET WINDUP INHIBIT	ok [lf]
UWF n [cr]	"n" is reset windup inhibit value as a percentage of prop set UNWIND FACTOR	portional band ok [lf]
ICP n [cr]	set INTEGRAL CLIP	ok [lf] I to cut
DER n [cr]		ok [lf]
IDB n [cr]		ok [lf]
MNO n [cr]		ok [If] ercentage of full power
MXO n [cr]		ok [If] <i>III power</i>
Use this command to RETURN the se	elected 600/620 Controller to FRONT PANEL oper	ATION:
LCL [a]	place selected unit under LOCAL control (Front	: Panel) ok [lf]

▲ MANUAL CONTROL

Use these commands to execute real-time control of setpoints and event outputs of your Model 600A/620A Environmental Chamber System Controller.

СОММАND	FUNCTION	R E P L Y
1 Use these commands to first select the AOO	7620 YOU WISH TO ADDRESS AND PLACE THE SELECTED	LINIT LINDER COMPLITER CONTROL
ADR [cr]	read the ADDRESS of the selected unit	
#nn [cr]	select the 600/620 ADDRESS	[no reply]
IFS [cr]	request INTERFACE STATUS	LCL [IF] or RMT [IF]
RMT SELECT [cr]	change to REMOTE mode (computer control)	ok [lf]
	To transfer control back to the 600/620 Front Panel, use the	command LCL [cr]
2 Use this command to place BOTH CHANNELS	OF THE SELECTED UNIT INTO THE MANUAL SETPOIN	T MODE:
CH1 M [cr]	set Channel 1 Setpoint operation to MANUAL Mode	ok [lf]
	Use the command "CH2 M" to place Channel 2 Setpoint opera	tion under Manual control
A LICE THESE COMMANIES TO DEAD THE STATUS OF S		
3 USE THESE COMMANDS TO READ THE STATUS OF E	read the MANUAL SETERING for Channel 1	vvv v [lf]
נון	read the MANUAL SETTOINT for Channel 2	xxx x [If]
PR1 [rr]	read the current PROCESS VALUE for Channel 1	xxx x [lf]
PR2 [rr]	read the current PROCESS VALUE for Channel 2	xxx x [lf]
	read the current EVENT STATUS	1•••5•7•A•••F••H [lf]
[0]	Events that are "ON" are noted by identifying character; events	s that are "OFF" are noted by "•"
4 Use these commands to EXECUTE REAL-TIME	CONTROL OVER SETPOINTS, EVENTS AND THERMOE	Boost of the selected unit:
CM1 [cr]	clear the MANUAL SETPOINT for Ch 1 (Use "CM2"	for Ch2) ok [If]
MS1 nnn.n [cr]	set the MANUAL SETPOINT for Ch 1 (Use "MS2"	for Ch2) ok [If]
CH1 % OUTPUT [cr]	set Ch1 to %OUTPUT (PCT) mode (Use "CH2" fo	r Ch2) ok [lf]
M%1 nnn [cr]	set nnn% output on Ch1 (Use "M%2" for Ch2) PCT OUTPUT mode must be selected	ok [lf]
MBD nn [cr]	set THERMOBOOST value for COOL action	ok [lf]
MBI nn [cr]	set THERMOBOOST value for HEAT action	ok [lf]
MA1 nn [cr]	set ALARM LIMIT for Ch 1 (Use "MA2" for Ch2) Limit defined is deviation from manual setpoint	ok [lf]
MEN n n [cr]	set specified events to "ON"	ok [lf]
MEC [cr]	clear all MANUAL events	ok [lf]
MEF n n [cr]	set specified events to "OFF"	ok [lf]
5 Use this command to RETURN the selected 6	00/620 Controller to FRONT PANEL OPERATION	
LCL [cr]	place selected unit under LOCAL control (front panel	l)ok [lf]

I E E E 4 8 8O P T I O N SRefer to the back cover of this guide for a listing of additional IEEE488 instructions.

▲ CONTROL PROGRAMMER

RUN, STOP (HOLD), RESET, ETC.

Use these commands to control operation of your 600A/620A Programmer – Run, Stop, Reset, etc. – from your computer or smart terminal.

C U M M A N D	FUNCTION RE	ΡΙΥ
1 Use these commands to first select the 600	1/620 you wish to address and place the selected unit linder comput	FR CONTROL:
ADR [cr]	read the ADDRESS of the selected unit	nn [lf]
#nn [cr]	Use this QUERY command if your installation has two or more 600/620s	[no ronly]
IFS [cr]	request INTERFACE STATUS	[f] or RMT [lf]
	Control is either from the Front Panel ("LCL") or from your computer ("RMT")	
RMI SELECI [cr]	To transfer control back to the 600/620 Front Panel, use the command LCL [cr]	ok [If]
2 Use these commands if you wish to READ th	ie selected program and step or SELECT a specific program and step:	
PRN [cr]	read the selected PROGRAM number	nn [lf]
SIN [Cʃ] DDN nn [cr]	solect a PPOCPAM number	nnn [lf] ok [lf]
STN nnn [cr]	select a STEP NUMBER in the selected program	ok [lf]
³ Use these commands to RUN a Program:	1.5	
RUN p [cr]	. run PROGRAM "p"	ok [lf]
RUN p s [cr]	. run PROGRAM "p" from STEP "s" "s" must be a RAMP/SOAK (STD) step	ok [lf]
4 Use these commands to STOP (HOLD) or R	ESET THE PROGRAM THAT IS CURRENTLY EXECUTING:	
STP [07]	STOP (HOLD) the program that is executing	ok [lf]
	. RESET THE STOPPED PROGRAM TO STEP 1	
VALUE)	RAMMED VALUES FOR THE CURRENT STEP: (NOTE: SEND COMMAND W7 O DA	IA IO READ CURRENI
, CS1 nnn[cr]		ok [lf]
CA1 nnn[cr]		ok [lf]
CCT hh mm ss[cr]	. Change Current TIME REMAINING in Step	ok [lf]
CCE 178AGH[cr]	. Change the Current EVENTS Note: all 16 events must be entered in proper sequence	ok [lf]
CCL nnnn[cr]	. Change Current LOOPS	ok [lf]
CCC nnnn[cr]	Change Current CYCLES	ok [lf]
	Change Current PROGRAMMED HEAI (Increase Thermal Boost)	OK [If]
		UK [II]
DIN [cr]	RESUME execution of program	ok [lf]
	PROGRAM # and STEP # must have been previously selected	טא נוון
7 Use this command to RETURN the selected	500/620 Controller to FRONT PANEL operation:	
LCL [cr]	place selected unit under LOCAL control (front panel)	ok [lf]

🔺 UPLOAD / DOWNLOAD

Save and Restore Programs and PID Values

Use these commands to transfer programs and tuning values between your computer and your Model 600A/620A Environmental Chamber System Controller.

C O M M A N D	FUNCIION	REPLY
Use these commands to first select the 600/	620 you wish to address and place the sel	ected unit under computer control:
ADR [cr]	read the ADDRESS of the selected unit	nn [lf] more 600/620s
#nn [cr]	select the 600/620 ADDRESS	[no reply]
IFS [cr]	request INTERFACE STATUS Control is either from the Front Panel ("LCL") or from you	LCL [IF] or RMT [IF]
RMT SELECT [cr]	change to REMOTE mode (computer control) To transfer control back to the 600/620 Front Panel, us	ok [If] the command LCL [cr]
2 Use these commands if you wish to SAVE (up	load) programs and pid parameters to the co	IMPUTER:
ULP n [cr]	save PROGRAM "n" from the selected unit The data for program "n" is represented by "xxxxxx.	×XXXXXXXXXX\\\ [lf]
ULD P [cr]	save ALL PID PARAMETERS from the selected u The PID data is represented by "xxxxxxxxxx"	ınit xxxxxxxxxx ok [lf]
ULD S [cr]	save the TOTAL PROGRAM MEMORY from the The data for all 200 steps of program memory is represe	selected unit . xxxxxxx [If] [If] ok [If] ented by "xxxxxxxxx"
3 Use these commands if you wish to RESTORE	(DOWNLOAD) PROGRAMS OR PID VALUES TO THE	selected 600/620 Controller:
DLP nn [cr] xxxxxxxx [cr]	Load PROGRAM "nn"	ok [If] viously saved using the "ULP nn" command
DLD P [cr] XXXXXXXX [cr]	Load PID Parameters	ok [If] using the "ULD P" command
DLD S [cr] XXXXXXXX [cr]	Load TOTAL PROGRAM MEMORY (all 200 step "xxxxxxxx" represents the total program memory pre	s)ok [If] viously saved using the "ULD S" command
4 Use this command if you wish to VERIFY a p	Rogram in the memory of the selected 600 /	'620 Controller:
VFP nn [cr] xxxxxxxx [cr]	VERIFY Program "nn" "xxxxxxxx" represents the data for Program "nn" prev Response is *ok [If] if verification is successful; otherwis	*ok [If] or *? [If] viously saved using the "ULP nn" command e response is *? [If]
5 Use this command to RETURN the selected 60	00/620 Controller to FRONT PANEL OPER	ATION:
LCL [cr]	place selected unit under LOCAL control (front	panel)ok [lf]

▲ IEEE-488 COMMAND SET

COMMAND

Use these commands to perform the specified IEEE-488 operations in connection with your Model 600A/620A Environmental Chamber System Controller.

FUNCTION

REPLY

RQN [cr]	ENABLE IEE488 service request	ok [lf]
RQD [cr]	DISABLE IEEE488 service request	ok [lf]
SRN n n [cr]	enable individual SRQ MASK BIT	ok [lf]
	The SRQ bits are assigned as follows:	
	Bit 0: Deviation Alarm Channel 1 or Channel 2	Bit 1: Event F (OFF to ON Transition)
	Bit 2: Setpoint Alarm Channel 1 or Channel 2	Bit 3: Input at TB-6, Pin #8
	Bit 4: Channel 1 Process Limit	Bit 5: Channel 2 Process Limit
	Bit 6: Serial Poll (Set Anytime SRQ Issued)	Bit 7: Unused
SRS [cr]	read IEEE488 SRQ mask status register	1••4••••[lf]
SRS •1•3•5•• [cr]	store image of the total SRQ mask	ok [lf]
SRC [cr]	clear enabled SRQ mask bits	ok [lf]
GTL [cr]	go to local	no reply
DCL n [cr]	clear device	no reply



From Computer PS/22	F	RS422 H	ookup		
Port, or JCS Model	Contro	ller #01	Control	ler #02	Controller #nn
404 Level Shifter	P6	J5	P6	J5	P6 J5
Shield (ground one end only) Send Receive Signal Ground (SG) Send + Receive +	1 1 1 1 1 1 1 1 1 1 1 1 1 1	entions (1) (1) (1) (1) (1) (1) (1) (1) (1) (1)	(1) Conne (4) Conne (6) (1) (2) (2) (2) (2) (2)	mal tettoms	Internal 1 Connections 1 4 - 6 - 9 - - - - - - - - - - - - - - - - - - - - - - - - - - - - - - -

VServices Inc. P.O. Box 157 Elk River, MN 55330

Phone: 763-241-1456 Fax: 763-241-1829 Email: jcsystems@tmcservices.net Web: www.jcsystemsinc.com